

# Capp *Rondo*

3 D SHAKER

**Capp *Rondo* 3D Shaker** operates at a variable speed range from 10 to 120 RPM. The brushless DC motor ensures stable, quiet and maintenance free operation, which is perfectly suited for extended runs. The three dimensional rocking and tumbling motion provides thorough homogenous mixing.

The user can save up to 9 individual protocols or use the pulse mode with Run-Pause-Run operation at a programmable time interval ranging from 30 to 99 seconds.

**Capp *Rondo* 3D Shaker** can be programmed to run from 1 up to 999 minutes, or continuously in the infinite mode.

To assure total stability on the bench, **Capp *Rondo* 3D Shaker** is equipped with non-slip rubber feet and comes with a spiked rubber mat. Optionally, a sticky pad mat can be ordered.

**Capp *Rondo* 3D Shaker** is easy to clean and complies with Protection Class IP 21.



## Specifications:

Weight (Capacity)	5kg
Speed Range	10-120 RPM
Motor Type	BLDC (Brushless DC motor)
Movement	Tumbling 7°
Run Time	1-999 minutes or infinite
Programmable Mode	Up to 9 user defined protocols
Pulse Mode	Run-Pause-Run programmable from 30 to 99 secs.
Display	7 Segment LED
Accessories	Spiked rubber mat (included), Sticky pad mat (optional)
Plate size	350x350mm
Protection Class	IP 21

## Ordering Information

Cat No.	Description
CR3-08X	CappRondo 3D Shaker 5kg capacity

[www.capp.dk](http://www.capp.dk) – +45 6613 6140 – [info@capp.dk](mailto:info@capp.dk)