



CAPP APS
 LANDBRUGSVEJ 10
 DK-5260 ODENSE
 DENMARK

TEL. (+45) 6613 6140
 FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK
 WEB: WWW.CAPP.DK



Capp *Rondo*

M A G N E T I C S T I R R E R

The Capp Rondo Magnetic Stirrer is a compact device that works at a maximum speed of 1500 rpm. It is manufactured with robust stainless steel top surface to resist moderate contamination of harsh chemicals. The adjustable speed and timer function can be set and monitored on a modern real-time digital display.

Not just a simple Stirrer:

- Adjust speed and timer.
- Last run memory feature.
- Pulse mode works for 30 seconds on each direction, for an improved mixing.
- Water spill resistant – IP 65 compliant.

Innovative design:

- Slim profile – only 16.5mm thick.
- Touch control pad and digital display.
- Flat surface for ease of cleaning.
- Maintenance free.

Important note: to maximize the performance of your Capp Rondo Magnetic Stirrer, Capp recommends that you use our original Teflon coated, 25mm high strength stirring bar. One unit is supplied with each instrument. Contact your local distributor to order more units.

Specifications	Capp Rondo Magnetic Stirrer
Maximum speed	1500 rpm
Speed setting	15 to 1500 rpm (increments of 50 rpm)
Load capacity	800 mL
Recommended stir bar length	25 mm
Dimensions	124 x 191 x 16.5mm
Weight	600 grams
Pulse function	YES
Timer function	1 to 99 minutes and infinite mode



Ordering Info

Cat No.	Description
CRS-15X	Capp Rondo Magnetic Stirrer 1500rpm w/ 1 stirring bar